

TRAUMA KIT NOW! ADVANCED KIT SUPPLY REFILL - TRAUMA KIT NOW! ADVANCED SUPPLIES REFILL

The Trauma Kit refill is not just another small tactical first aid kit. Each component was specifically chosen for its combat proven effectiveness combined with never before seen packaging techniques to make them as small as possible. The most popular MTKN configuration in service with police, military and medical providers worldwide. The ADVANCED kit is designed to treat traumatic injuries as per the MARCH protocol for those providers trained in TCCC concepts. The lightest and smallest fulltreatment trauma kit available. The Trauma Kit supplies were designed from the ground up around a specific need for an ultra-concealable trauma kit. Refill Kit Includes (1x Pair) Trauma Gloves (1X) Quikclot Combat Gauze (1X) Emergency Trauma Dressing 6" (1X) Hypothermia Blanket - Green / Silver (1X) Flat Fold Tape (1X) Compressed Gauze (1X) Trauma Shears (1X Twin Pack) Full Size Chest Seal (1X) Nasopharyngeal Airway (1X) Decompression Needle



Attributes

- Name: TRAUMA KIT NOW! ADVANCED SUPPLIES REFILL
- Manufacturer: BLUE FORCE GEAR
- Product no.: 430100763
- Mfr. No.: TKR-MED-ADV
- Size: 1 Person
- Delivery weight: 0.204kg
- Shipping height: 12mm
- Shipping width: 191mm
- Shipping length: 235mm
- UPC: 810073650848

Item details

Made in USA

Table of Contents

- [Startpage](#)
- [English: Trauma Kit Now! Advanced Kit Supply Refill Safety Instructions](#)
- [Italiano: Guida alle Istruzioni di Sicurezza per il Trauma Kit Now! Advanced Kit Supply Refill](#)

Trauma Kit Now! Advanced Kit Supply Refill

Safety Instructions

Introduction

Thank you for choosing the Trauma Kit Now! Advanced Kit Supply Refill. This guide provides essential safety instructions to ensure the safe and effective use of the kit. Please read this manual thoroughly before use to understand the proper precautions and procedures.

General Safety Guidelines

- Ensure the kit is stored in a cool, dry place away from direct sunlight.
- Regularly check the contents of the kit to ensure all items are in good condition and not expired.
- Keep the kit out of reach of children unless under adult supervision.
- Report any unsafe products or incidents to the relevant authorities.
- Stay informed about product recalls via the EU's Safety Gate platform.

Specific Safety Precautions for Use

- Always wear protective gloves when handling blood or bodily fluids.
- Use trauma shears with caution to avoid injury.
- Follow the MARCH protocol for treating traumatic injuries:
 - M Massive Hemorrhage
 - A Airway
 - R Respiratory
 - C Circulation
 - H Hypothermia
- Use the decompression needle only if trained and certified in its use.
- Dispose of used materials, such as gauze and gloves, in accordance with local regulations.

Instructions for Installation and Usage

1. Preparation:

- Ensure you are in a safe environment before beginning any first aid procedures.
- Gather all necessary items from the kit.

2. Using Trauma Gloves:

- Open the glove package and put on the gloves before handling any injured person.
- Dispose of gloves properly after use.

3. Applying Quikclot Combat Gauze:

- Remove the gauze from its packaging.
- Apply direct pressure to the wound using the gauze until bleeding stops.

4. Using the Emergency Trauma Dressing:

- Open the dressing package.
- Place the dressing over the wound and secure it in place.

5. Using the Hypothermia Blanket:

- Remove the blanket from the package.
- Wrap the injured person to maintain body temperature.

6. Using Flat Fold Tape and Compressed Gauze:

- Use the flat fold tape to secure dressings in place.
- Apply compressed gauze to wounds as needed.

7. Using Trauma Shears:

- Carefully cut clothing or other materials away from the wound site.
- Ensure the shears are stored safely after use.

8. Using the Full Size Chest Seal:

- Place the chest seal over the wound, ensuring it is airtight.
- Monitor for signs of respiratory distress.

9. Using the Nasopharyngeal Airway:

- Only use if trained; insert gently into the nostril until it reaches the back of the throat.

10. Using the Decompression Needle:

- Only use if trained; insert into the chest cavity as needed for tension pneumothorax.

Disposal Instructions

- Dispose of used gloves, gauze, and other single-use items in accordance with local regulations.
- Follow local guidelines for the disposal of medical waste.
- Keep the kit clean and restock it as necessary to maintain readiness.

Contact Information for Further Support

For any safety inquiries or further support, please reach out to the designated contact point within the EU. Ensure you have your product details ready for assistance.

Thank you for prioritizing safety with the Trauma Kit Now! Advanced Kit Supply Refill. Your preparedness can make a significant difference in emergency situations.

Guida alle Istruzioni di Sicurezza per il Trauma Kit Now! Advanced Kit Supply Refill

Introduzione

Benvenuto nella guida alle istruzioni di sicurezza per il Trauma Kit Now! Advanced Kit Supply Refill. Questo kit è progettato per fornire forniture di pronto soccorso avanzate per il trattamento delle lesioni traumatiche. È fondamentale seguire le linee guida di sicurezza per garantire un uso corretto e sicuro del prodotto.

Linee Guida Generali di Sicurezza

- Assicurati che il kit sia sempre conservato in un luogo fresco e asciutto.
- Controlla regolarmente il contenuto del kit per assicurarti che le forniture siano complete e non scadute.
- Usa il kit solo per le indicazioni previste e secondo le istruzioni fornite.
- Non utilizzare forniture danneggiate o compromesse.
- In caso di emergenza, contatta i servizi di emergenza locali oltre a utilizzare il kit.

Precauzioni di Sicurezza Specifiche per l'Uso

- Indossa sempre i guanti per traumi forniti quando si maneggiano ferite o sangue.
- Utilizza la Quikclot Combat Gauze per il controllo delle emorragie, seguendo le istruzioni sulla confezione.
- Assicurati che la benda traumatica d'emergenza sia applicata correttamente per evitare ulteriori danni.
- Utilizza la coperta ipotermica per mantenere la temperatura corporea della persona in caso di shock.
- Le forbici per traumi devono essere utilizzate con attenzione per evitare infortuni accidentali.

Istruzioni per l'Installazione e l'Uso

- Apri il kit e verifica la presenza di tutte le forniture elencate.
- In caso di lesioni traumatiche, segui il protocollo MARCH:
 - M: Controllo delle emorragie.
 - A: Apertura delle vie aeree.
 - R: Respirazione.
 - C: Circolazione.
 - H: Ipotermia.
- Inizia con il controllo delle emorragie utilizzando la Quikclot Combat Gauze e la benda traumatica.
- Se necessario, utilizza la via aerea nasofaringea per garantire una corretta respirazione.
- Dopo aver fornito assistenza, monitora costantemente la persona fino all'arrivo dei soccorsi.

Istruzioni per lo Smaltimento

- Smaltisci le forniture scadute o danneggiate in conformità con le normative locali sullo smaltimento dei rifiuti.
- Non gettare le forniture nel normale cestino se contengono materiali biologici.
- Considera di contattare un servizio di smaltimento speciale per i materiali pericolosi, se necessario.

Informazioni di Contatto per Ulteriore Supporto

Per domande o dubbi relativi alla sicurezza del prodotto, ti invitiamo a contattare il tuo rivenditore o il produttore. Assicurati di avere a disposizione il numero di lotto e la data di scadenza del prodotto per facilitare il supporto.

Conclusione

Seguire queste linee guida di sicurezza è essenziale per garantire un uso efficace e sicuro del Trauma Kit Now! Advanced Kit Supply Refill. In caso di dubbi o domande, non esitare a cercare assistenza. La tua sicurezza è la nostra priorità.